

LabJack

Published on *LabJack* (<https://labjack.com>)

[Home](#) > [Support](#) > Modbus TCP Example Code Table

Modbus TCP Example Code Table

Modbus TCP Example Code Table (Referencable)

Direct Modbus TCP Example Code



T7, T4

<u>Libraries...</u>	<u>Windows</u>	<u>Mac</u>	<u>Linux</u>
<u>LabVIEW</u>	✓	✓	✓
<u>Node.js</u>	✓	✓	✓
<u>Python</u>	✓	✓	✓
<u>C/C++ (Streaming)</u>	✓	✓	✓

To use a 3rd party Modbus program, such as a Rockwell/Allen Bradley application, follow the instructions on the [Modbus Client Applications](#) page.

More information about Modbus as well as some getting started information can be found on the [Modbus API Documentation](#) page. A list of the available Modbus registers is available on the [Modbus Map](#) page.
