

# LabJack

Published on *LabJack* (<https://labjack.com>)

[Home](#) > [Products](#) > U6

---

## U6

## U6

USB multifunction DAQ device with 14 analog inputs (16- to 18- bit), 2 analog outputs (12-bit), 20 digital I/O, 2 counters, and 4 timers.

U6

### Quantity

\$345.00

Quantity Discounts:

1 - 9 = \$345.00 per unit

10 - 29 = \$317.40 per unit

30 - 99 = \$296.70 per unit

100 - 499 = \$276.00 per unit

500+ = \$262.20 per unit

 [U6 Product Brief](#)

**Stock Status:** In Stock

## USB Multifunction DAQ

The U6 family is great for acquiring signals from thermocouples, strain gages, bridge circuits, encoders, and much more. With instrumentation amplified inputs, and software selectable gain settings, the U6 provides a cost effective method to measure analog signals which can't be resolved by traditional 12-bit DAQ products. Common applications include laboratory research, industrial control and monitoring, and prototype development.

▶ [Product Highlights](#)

▶ [Hardware Overview](#)

▶ [Product Variations](#)

▶ Other Highlights

▶ Package Contents

▶ Accessories

## Support Resources

▶ Getting Started

▶ Software Options

▶ Datasheet and Application Notes

▶ Contact Us

---